

Enhancing Upstream Processes for High-Yield, High-Quality AAV Vector Production Using a Novel In-House Cell Line

ASGCT 2025

Dicky Gilmore

Sr Scientist

Upstream Process Development

OXB; Bedford, MA



Legal disclaimer

This presentation does not constitute an offer to sell or a solicitation of offers to buy Ordinary Shares (the “Securities”). Although reasonable care has been taken to ensure that the facts stated in this presentation are accurate and that the opinions expressed are fair and reasonable, the contents of this presentation have not been formally verified by Oxford Biomedica plc (“OXB” or the “Company”) or any other person. Accordingly, no representation or warranty, expressed or implied, is made as to the fairness, accuracy, completeness or correctness of the information and opinions contained in this presentation, and no reliance should be placed on such information or opinions. Further, the information in this presentation is not complete and may be changed. Neither the Company nor any of its respective members, directors, officers or employees nor any other person accepts any liability whatsoever for any loss howsoever arising from any use of such information or opinions or otherwise arising in connection with this presentation.

This presentation may contain forward-looking statements that reflect the Company's current expectations regarding future events, its liquidity and results of operations and its future working capital requirements. Forward-looking statements involve risks and uncertainties. Actual events could differ materially from those projected herein and may depend on a number of factors.



Presentation Overview

1

OXB's AAV Platform

2

**Approach to Upstream
Process Enhancement**

3

**Developing an In-House
Cell Line for AAV**

4

**'Plug-n-Playability' of
new Upstream Platform**





Unique competitive positioning

- ✓ **Best-in-class capabilities** across AAV, lentivirus & other vector types
- ✓ **Trusted by global industry leaders** – successful collaborations with big pharma, established biotech and emerging biotech
- ✓ **State-of-the-art facilities & scalable production capabilities** designed to meet the growing demand for C>s
- ✓ **Deep scientific know-how** – a team of world-leading specialists in viral vector optimisation
- ✓ **Cutting-edge technology** – leveraging 30 years of insights to enhance speed, efficacy, quality and safety in new therapies
- ✓ **Global reach & strategic positioning** with manufacturing facilities located in key biotech hubs



30
Years of
manufacturing
experience



960+
Successful
GMP batches
since 2014




45+
Client
programmes




30+
IND
submissions



65+
Successful
audits



OXB's AAV Platform



Upstream Optimization Focuses

1

Cell Line

- Current platform uses commercially available HEK293 cell line for highly robust, high-titer AAV production
- **Goal: develop in-house cell line with equivalent productivity and product quality**




2

Transfection and Production

- Always evaluating new materials and methods to remain best-in-class
- Transfection practices and parameters are robust, cost-effective, and scalable
- **Plug and Play!** Platform must be optimized to fit both triple and dual transfection across all capsid serotypes

Creating an Enhanced, Next-Generation Upstream Process

Conducted process development work on multiple Upstream factors, notably:



Transfection Reagent

Bench-scale screening to identify best commercially available reagent

Top candidate selected for **titer, percent full capsids, reduced batch costs, and transfection complex stability**




Transfection Parameters

Full factorial DOE to optimize:

- **Transfection reagent : DNA ratio**
- **DNA : Cell ratio**
- **Transfection Cell Density**

Ensured setpoints and ranges appropriate for large scale MFG



Process Additives


Identified commercially available additive to boost vg productivity

Optimized **addition amount** and **addition timing** for maximum yield and product quality


3.2x Titer Improvement from Upstream Process Optimization

AAV9 Dual Transfection; 2L Bioreactor Scale

Titer



Comparable packaging seen with each process improvement




Plug-and-Playability


Will this new Upstream platform fit multiple cell lines and serotypes?

3.2x Titer Improvement from Upstream Process Optimization

AAV9 Dual Transfection; 2L Bioreactor Scale



Comparable packaging seen with each process improvement




Plug-and-Playability


Will this new Upstream platform fit multiple cell lines and serotypes?

Creation of an Internal OXB Cell Line for AAV Production

Parental HEK293T Cell Line



Used in OXB's
Lentivirus platform



SV40T
Removal
and
Screening

Final OXB Cell Line





Whole Genome Sequencing

No SV40T or
Mad7 Detection


Enhanced Upstream Process Improves Titer in Both Cell Lines

AAV9 Dual Transfection; 2L Bioreactor Scale



Enhanced Upstream Process Improves Titer In Multiple Serotypes

Increased performance from dual transfection seen in **both cell lines and platforms**




Process C: External Cell Line
Process D: Internal Cell Line
All 2L Bioreactor Scale Data

2e15 vg/L or higher in all 6 tested serotypes for dual transfection

Good Product Quality Maintained Regardless of Process or Cell

Final DS		Process C		Process D
		External Cell Line	HEX2s52	HEX2s52
VP Ratio		1:1:10	1:1:8	1:1:8
AUC	%Empty	13.5%	7.2%	10.0%
	%Partial	7.2%	12.1%	10.0%
	%Full	79.3%	80.7%	80.0%
Aggregation		0.8%	0.2%	0.4%
Host-Cell Protein (ng/1E13 vg)		BLOQ	BLOQ	BLOQ
hcDNA (ng/1E13 vg)		6.4	25.7	4.6
Residual RepCap DNA (copies/1E13 vg)		4.3E9	BLOQ	5.2E9
Residual pHelper DNA (copies/1E13 vg)		5.0E8	4.5E8	2.1E9
Residual Ad E1A DNA (copies/1E13 vg)		BLOQ	BLOQ	BLOQ

Potency comparable or improved in all serotypes for Process D



Summary

1

Upstream platform enhanced via new transfection materials and parameter optimization

2

In-house cell line engineered to fit upstream process and produce equivalent titer and PQ

3

Optimized Upstream “Process D” shown to be plug-and-play with multiple cell lines and serotypes, consistently achieving $>2e15$ vg/L



Acknowledgement

Thank you

Upstream Development

- Shannon Frost
- Marissa Stanvick
- Diane Golebiowski
- Chase Rebidue
- Laura Van Lieshout

Cell Line Development

- Hannah Stewart
- Louis Frost
- Magda Reis
- Rebecca Astley

Analytical Development

- Serena Dollive
- Jeff Gagnon
- Kiran Adhikari
- Eugenia Fandunyan

Downstream Development

- Alex Meola
- Carlos Chong
- Nishadi Gamage
- Alejo Zacarias
- Danny Lassiter
- Iraj Ghazi
- Vincent Thari

Additional Thanks

- Ify Iwuchukwu
- Jin Yin
- Scott Barbick
- John Maravich
- Kyriacos Mitrophanous

Let's deliver life-changing therapies together

A global quality and innovation-led
CDMO in cell and gene therapy

[OXB.com](https://oxb.com)

